



INSTRUCTIONS FOR USE

BUNION PROTECTOR

INDICATION

Using the GO Medical® bunion protector relieves the pain associated with the pressures and frictions on the bunion.

CONTRA-INDICATIONS

Do not use on injured skin. In case of arteritis, diabetes, lower limb neuropathy or ulcers, seek advice from your doctor or podiatrist before use and then regularly monitor your foot.

PRECAUTIONS FOR USE

- The hallux valgus protection pad should be worn next to the skin, on the right or left foot.
- Slip your foot into the protection pad, putting your great toe in the tubular part. The gel module should be placed on the hallux valgus.
- To help the protection pad to remain in position, wear tights, socks or stockings over it.

FEATURES

For bunion protection, size matters: carrying a protection that is too bulky in the shoe can aggravate the problem further or create undesirable pain on other areas of the foot. The GO Medical® protector is therefore thin and lightweight, but still very effective. The silicone gel helps avoid frictions

on the bunion and was deliberately limited to 1mm of thickness.

DESIGN

Thin and lightweight, the protector can be worn in any shoe.

ADVICE

If your great toe is still mobile laterally, we recommend wearing the GO Medical® flexible bunion corrector. It limits the evolution of the deformation, relieves joint pain and avoids rubbing/friction on the bunion.

CLEANING & LIFESPAN

Wash in a machine at 30°C / 86°F or by hand with soap and water (without rubbing) and leave to dry naturally. Do not iron. Reusable for several months.

COMPOSITION: 61% polyamide/elastane fabric, 39% silicone.



Machine washable at 30°C / 86°F.



Do not iron or steam.



Do not dry clean.



Do not bleach.



Do not tumble-dry.



See packaging.



Made in China
Imported By ProBrace



www.probrace.com