Ref.: AB100

Boxia® Plus

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Ref.: AB100

Boxia® Plus

INSTRUCTIONS FOR USE, STORAGE AND WARRANTY

Dear Customer,

Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If you have any questions or concerns, please contact your doctor, orthopaedic specialist or our customer service department.

ORLIMAN S.L.U. guarantees all its products as long as the original configuration has not been manipulated or altered except for the intended use as described in these instructions.

If the products are used in combination with other products, replacement parts or systems, make sure they are compatible and made by Orliman. It does not guarantee any products with altered characteristics due to improper use, defects or breakage of any kind. The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country.

Orliman would like to thank you for choosing this product and hopes you a speedy recovery.

REGULATIONS

MD This article is defined as a class I medical device. A Risk Analysis (UNE EN ISO 14971) has been carried out, minimising the existing risks. Tests have been in accordance with European Regulation UNE-EN ISO 22523 on Prostheses and Orthoses.

INDICATIONS

- Weakness or inability to dorsiflex the ankle (dropfoot) as a result of neurological injuries, congenital alterations or muscular disorders.
- Flaccid paralysis associated with hemiparesis.

FITTING INSTRUCTIONS

For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended. If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

When fitting the product, you must adhere to the following instructions:

Using an interface (i.e. a sock) is recommended.

A ADJUSTING THE PERIMETER

- The product can be adapted to a specific ankle perimeter by adjusting the part containing the elastic guide (pin).
- Depending on the patient's ankle perimeter in the malleolus area, position this part between the first and last of the 4 holes on the supramalleolar part. See the following table for guidance:

Size	(cm)	(cm)	(cm)	(cm)
1	14	16	18	20
2	20	22	24	26
3	26	28	30	32

B SECURING THE HOOKS

- If using the anti-slip part with a hook, until the laces on the footwear and position it between the shoe tongue and laces so it is centred. Secure the lower micro-hook by holding the laces.
- If using the "S" hooks, you must insert them in the eyelets with the longest part inside the shoe. The elastic traction strap must go through both hooks.

C ADAPTING THE SUPRAMALLEOLAR PART

- Open the supramalleolar part until you hear a click; the part will be flat.
- Position the supramalleolar part in the back of the ankle so it coincides with the centre of the silicone pad against the Achilles tendon.
- Press the ends of the part toward the leg until you hear a click which unlocks the part to be wrapped around the ankle.
- Secure the tip of the widest micro-hook so the buckle is between the 2 narrowest straps.
- Pull and secure the narrowest straps.
- Insert the elastic traction band through the guide (pin) on the front.
- Secure the elastic band to the hook (or hooks) and pull on the ends with the micro-hook and then secure them to get the desired traction.

D FITTING THE INNER STRIP Ref.: AB12/AB22

- Secure the strip in the middle of the foot with the Velcro closures in the instep. Make sure the closure hook is in the right direction. Then, insert the elastic traction band through the guide on the front of the supramalleolar part. Secure it to the hook and pull on the ends to get the desired traction.

E INSERTING THE SHIN SUPPORT Ref.AB14/AB24

- The shin support must correspond to the size of the Boxia AB100, after being previously adapted and checked for the patient.
- Secure using the micro-hook on the inside and front of the supramalleolar part.
- The "Y"-shaped shin support has 3 ends. The shortest one along with the longest and narrowest will hug the shin at the height of the top of the fibula. The remaining end will go over the front of the tibia to be secured to the supramalleolar part on the inside.


PRECAUTIONS

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.

To avoid minor discomfort caused by sweating, we recommend using some type of cotton fabric to separate the skin from contact with the product material. For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

Products marked with the  symbol contain natural rubber latex and can cause allergic reactions in people sensitive to latex.

Products marked with the  symbol contain ferromagnetic components and, therefore, extreme precaution must be taken if you undergo an MRI scan or are exposed to radiation associated with diagnostic or therapeutic procedures.


RECOMMENDATIONS-WARNINGS

The use of these products is conditioned by the indications. Although the product is not defined as a single-use device, using it on a single patient only is recommended and only for the intended purposes as described in these instructions or by a healthcare professional.

When disposing of the product and its packaging, you must strictly adhere to the legal regulations in your community.

RECOMMENDATIONS FOR STORAGE AND WASHING

When not using the product, store it in the original packaging in a dry place at room temperature. Stick the Velcro to each other (if the orthotic device has them), frequently wash by hand with warm water (30° C max.) and mild soap. To dry the product, use a dry towel to absorb as much moisture as possible and let it dry at room temperature. Do not hang it up or iron the product and do not expose it to direct heat sources such as stoves, dryers, direct sun exposure, etc. When using or cleaning the product, do not use abrasive or corrosive substances, alcohol, ointments or liquid solvents. If not dried off properly, the detergent residue may irritate the skin and cause the product to deteriorate.

EN The manufacturing date is included in the batch number which can be found on the packaging label as , in the following way: the second and third digits represent the year of manufacture and the fourth and fifth digits represent the month.

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