
ENGLISH

	Medical Device
	Caution: This device contains natural rubber latex which may cause allergic reactions.

INTENDED USE

The device is intended to provide gross immobilization of the lumbar spine and/or compression

Indications for use

For lower back conditions that may benefit from compression and gross motion restriction, such as:

- Lower back pain (Lumbago)
- Lumbosacral sprain/strain
- Degenerative joint disease

Contraindications

Must not be used by individuals for whom compression is contraindicated.

Warnings and Cautions:

The device should not be used for extended periods of time and removed during long periods of rest, e.g., sleep. The device should be used directly on intact skin and never over an open wound. Inspect the skin under the device's contact area before and after each use.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference. Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stay with the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device. The device is for single patient – multiple use.

USAGE

Size Selection

Measure the waist circumference. The sizing chart is shown on the package. If between sizes always select the smaller size for optimal fit.

Device Application

- Attach the SmartPad to the center of the loop surface of the device. The narrow side of the pad should face the lower edge of the device and the silicone dots should face the body. Press the device on to the fingers through the loops attached to the front ends of the belt (**Fig. 1**). The pad must sit in the middle of the back.
- Put the two ends of the belt forward evenly and guide them one over the other: press the left side against the stomach (**Fig. 2**) and attach the right side over the left (**Fig. 3**). When fastening the device, remove your fingers from the left loop.
- The device should fit snugly with high compression (**Fig. 4**). The SmartPad can be removed, if desired. When fitting with high compression, make sure that it is comfortable to wear.

Device removal

Pass your fingers through the loops and unfasten the belt ends from each other to remove the device from the body.

Cleaning and care

Wash regularly to retain the compression of the textiles.

- Machine-wash with a mild detergent at cold temperature (30°C)
- Air dry

Note: Secure the hook and loop fasteners before washing to avoid damaging the textiles.

Note: Do not tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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