ENGLISH

SYMBOLS

MD Medical Device

Caution: This device contains natural rubber latex which may cause allergic reactions.

INTENDED USE

The device is intended for compression and pain relief of the elbow

Indications for use

- Medial and lateral epicondylitis
- Elbow pain and swelling

Contraindications

Must not be used by individuals for whom compression is contraindicated.

- Warnings and Cautions:**
- The device should not be used for extended periods of time and removed during long periods of rest, e.g., sleep.
 - Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease, need to be extra careful and are advised to consult with a healthcare professional before using the device. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this device, discontinue the use immediately and consult with a healthcare professional.

GENERAL SAFETY INSTRUCTIONS

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

Read these instructions carefully before use. Keep them for future reference.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.
- The device is for single patient – multiple use.

USAGE

Size Selection

Measure the circumference 12 cm (4 3/4 in) below the elbow.

The sizing chart is shown on the package. If between sizes always select the smaller size for optimal fit.

Device Application

- Remove the RadialStrap from the device (**Fig. 1**).
- Fold the top edge of the device outward so that the FoldLock edge becomes visible (**Fig. 2**).
- Grasp the device near the SmartPad pocket and pull it over the elbow. Make sure that the outer and inner bony protrusions of the elbow (epicondyls) are positioned in the cut of the SmartPad (**Fig. 3**).
- Unfold the top part of the device so that the FoldLock edge is in contact with the skin (**Fig. 4**).
- Fasten the RadialStrap onto the device and wrap it from the inner side around the arm. Adjust the pressure intensity according to your needs (**Fig. 5**).

Device removal

Undo the RadialStrap and fold the top edge of the device outward so that the FoldLock edge becomes visible. Grab the device at the bottom edge near the SmartPad pocket and pull the arm.

Cleaning and care

Wash regularly to retain the compression of the textiles.

- Machine-wash with a mild detergent at cold temperature (30°C)
- Air dry.

Note: Secure the hook and loop fasteners before washing to avoid damaging the textiles.

Note: Do not tumble dry, iron, bleach, or wash with fabric softener. **Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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