



**PUSH YOUR
PERFORMANCE**

E⁺motion[®]

medi Performance Collection



medi. I feel better.

A person in athletic wear is shown in motion, running or jumping, against a background of a 3D grid pattern. The person is on the right side of the frame, with their right leg extended forward and their left leg bent. The grid pattern is composed of dark grey lines on a black background, creating a sense of depth and movement.

PUSH YOUR

PERFORMANCE

*WITH FUNCTION THAT HELPS YOU
ACHIEVE YOUR GOALS.*



E+motion®

For every level: E+motion offers personalised functional benefits – for different sporting goals and physical requirements. These supports unlock potential and motivate you to: Push your Performance.

WITH E+MOTION® PERFORMANCE IS CAPITALISED

The new generation features premium bandages for five joints and combines more than 70 years of compression experience with innovative technology made in Germany.

The triple combination of performance material and performance fit, as well as performance compression knit provides effective support for every movement. The symbiosis of a modern, sporty look and medical benefits makes the supports particularly suitable for active people who want more: more stability, comfort and performance.

Give your customers and patients function that helps – and benefit from the impressive potential of E+motion supports.



**THESE
3 FUNCTIONAL
FEATURES ARE WHAT
MAKE THE DIFFERENCE
AND CAN BE FOUND IN
ALL E+MOTION SUPPORTS**

1

PERFORMANCE MATERIAL

With merino wool and Dry Skin Technology.



MERINO

Is natural, odour-resistant and thermoregulating. As it is finished on the inside, there is no pilling effect – even after repeated washing.

MULESING-FREE

Mulesing-free wool eliminates the painful “mulesing” process for sheep – for greater animal welfare.



DRY SKIN TECHNOLOGY

Ensures effective moisture wicking. The support absorbs moisture up to twice as fast* and transports it to the surface where it evaporates.

2

PERFORMANCE FIT

Ergonomic sports knit – pleasant with every move.



ERGONOMIC SHAPE

For a virtually crease-free, secure fit.




XL MOVEMENT ZONES

Provide maximum flexibility in the flexion areas – for greater freedom of movement exactly where it is needed.

3

PERFORMANCE COMPRESSION

Stability, massage and power in the knit.



NOTICEABLY MORE STABILITY AND POWER IN THE KNITTED FABRIC

Highly elastic yarns with strong elastic recovery (strong natural tendency to return to their original shape).

The result: more than 25% tighter and 10% thinner than other knits* and a very compact knit pattern.



3D SILICONE QUARTZ PADS

Create a soothing massage effect during movement and stimulate circulation. Swelling can be relieved more rapidly.

*medi product comparison

Genumedi® E+motion®



Genumedi® E+motion®



PERFORMANCE MATERIAL

Merino wool: odour inhibiting and thermoregulating, mulesing-free.

Dry Skin Technology for optimal moisture management.

PERFORMANCE FIT

Perfect-Fit top band – fits with every move.

Grip-Ons for easy donning.

XL-movement zone for even better performance.

Ergonomic shape.

Perfect fit silicone pad for each size.

PERFORMANCE COMPRESSION

3D-silicone-quartz-pad for massage and reduction of swellings.

Perceptibly more stability and power* in the knit – for an even better body feeling.

Available in sizes 1 – 7

Aid number 05.04.01.0009

* Knitted fabric recovers original shape well



Levamed[®] F+ on[®]

Levamed® E+motion®

Supplied with
an application
aid



PERFORMANCE FIT

- Comfortable **sock cuff**.
- **Extra large stretch zone** in the instep area.
- Ergonomic shape.
- **XL movement zone** for an even better performance.

PERFORMANCE COMPRESSION

○ **Extra thin 3D silicone-quartz pad** for massage and reduction of swellings.

○ **Perceptibly more stability and power*** in the knit – for an even better body feeling.

PERFORMANCE MATERIAL

○ **Merino wool:** odour inhibiting and thermoregulating, mulesing-free.

○ **Dry Skin Technology** for optimal moisture management.

Available in sizes 1 – 5


Aid number 05.02.01.0002

* Knitted fabric recovers original shape well



Epicomed[®] otion[®]

Epicomed® E+motion®



PERFORMANCE FIT

- **Stretch ribbed cuff** – comfortable for each upper arm circumference.
- **Anti Slip System** – holds with every movement.
- **XL movement zone** for even better performance.
- **Ergonomic shape.**

PERFORMANCE COMPRESSION

○ **3D silicone-quartz pad** for massage and reduction of swellings.

○ **Removable band strap** with pressure scale.

○ **Perceptibly more stability and power*** in the knit – for an even better body feeling.

PERFORMANCE MATERIAL

○ **Merino wool:** odour inhibiting and thermoregulating, mulesing-free.

○ **Dry Skin Technology** for optimal moisture management.

Available in sizes 2 – 6

Aid number 05.08.01.1001

* Knitted fabric recovers original shape well

Manum d active E+moti

A woman with dark hair tied back, wearing a black sports bra and leggings, is captured in a dynamic pose while exercising with a thick, black battle rope. She is wearing a black wristband on her right wrist. The background is dark with a pattern of horizontal lines that create a sense of depth and movement. The lighting highlights her muscles, suggesting a high-intensity workout.

Manumed active E+motion®

PERFORMANCE COMPRESSION

Quartz strap for stabilization.

Relief for the ulnar styloid thanks to strap opening.

Perceptibly more stability and power*
in the knit – for an even better body feeling.

PERFORMANCE FIT

Ergonomic palm hand splint.

XL movement zone for
increased performance.

Seamless knit.

Ergonomic shape.

Supplied with
an application
aid



PERFORMANCE MATERIAL

Merino wool: odour inhibiting and
thermoregulating, mulesing-free.

Dry Skin Technology for optimal
moisture management.

Available in sizes 1 – 4

Aid number 05.07.02.0127

* Knitted fabric recovers original shape well



Lumbar[®] plus E+motion

Lumbamed® plus E+motion®



PERFORMANCE COMPRESSION

- **Breathable 3D-quartz-pad** for massage and decrease of muscle tension.

- **Perceptibly more stability and power*** in the knit – for an even better body feeling.

PERFORMANCE MATERIAL

Merino wool: odour inhibiting and thermoregulating, mulesing-free.

Dry Skin Technology for optimal moisture management.

PERFORMANCE FIT

- **XL movement zone** due to extra-large leg cut.

- **Performance fastenings** with doubled velcro and ergonomic hand loops.


- **Ergonomic shape.**

Unisex sizes 1 – 5
for efficient stock management
Aid number 23.14.03.2018

* Knitted fabric recovers original shape well




Genumedi® E+motion®



Ausführung / model	Item no.
Standard / standard	K.100.60*
Extraweit / extrawide (x-wide)	K.100.61*

Messung
seitlich innen
Measurement
on the medial
side



(d) cm	(f) cm	(f) cm x-wide	(d) inches	(f) inches	(f) inches x-wide	Größe / size
25 – 28	37 – 40	–	9 ¾ – 11	14 ½ – 15 ¾	–	*1
28 – 31	40 – 43	–	11 – 12 ¼	15 ¾ – 17	–	*2
31 – 34	43 – 46	46 – 49	12 ¼ – 13 ½	17 – 18	18 – 19 ½	*3
34 – 37	46 – 49	49 – 52	13 ½ – 14 ½	18 – 19 ½	19 ½ – 20 ½	*4
37 – 40	49 – 52	52 – 55	14 ½ – 15 ¾	19 ½ – 20 ½	20 ½ – 21 ¾	*5
40 – 43	52 – 55	55 – 58	15 ¾ – 17	20 ½ – 21 ¾	21 ¾ – 22 ¾	*6
43 – 46	55 – 58	58 – 61	17 – 18	21 ¾ – 22 ¾	22 ¾ – 24	*7


WHY CHOOSE GENUMEDI E+MOTION?

Breathable sports knit and perfect fitting pad with open edging to minimise retro-patellar pressure.

Aid number 05.04.01.0009



Levamed® E+motion®



Item no. K.010.60*

Messung
Measurement



cm	inches	Größe / size
17 – 20	6 ¾ – 8	*1
20 – 22	8 – 8 ¾	*2
22 – 24	8 ¾ – 9 ½	*3
24 – 26	9 ½ – 10 ¼	*4
26 – 28	10 ¼ – 11	*5


WHY CHOOSE LEVAMED E+MOTION?

Breathable sports knit and extra thin pad for (sports) shoes.

Aid number 05.02.01.0002




Epicomed® E+motion®



Item no. K.400.60*

Messung
Measurement



cm	inches	Größe / size
20 – 23	8 – 9	*2
23 – 26	9 – 10 ¼	*3
26 – 29	10 ¼ – 11 ½	*4
29 – 32	11 ½ – 12 ½	*5
32 – 35	12 ½ – 13 ¾	*6


WHY CHOOSE EPICOMED E+MOTION?

Breathable sports knit and individual pressure regulation with additional webbing strap with scale.

Aid number 05.08.01.1001



Manumed active® E+motion®



links / left	rechts / right
K.412.60*	K.413.60*


Messung Measurement



Handgelenk / wrist		
cm	inches	Größe / size
13 – 15	5 ¼ – 6	*1
15 – 17	6 – 6 ¾	*2
17 – 19	6 ¾ – 7 ½	*3
19 – 21	7 ½ – 8 ¼	*4




Lumbamed® plus E+motion®



Item no. K.210.60*


Unisex

Messung Damen Measurement ladies



Hüfte / hip		
cm	inches	Größe / size
76 – 86	30 – 34	*1
86 – 98	34 – 38 ¾	*2
98 – 111	38 ¾ – 43 ¾	*3
111 – 125	43 ¾ – 49 ¼	*4
125 – 140	49 ¼ – 55	*5

Messung Herren Measurement men



Taille / waist		
cm	inches	Größe / size
66 – 76	26 – 30	*1
76 – 88	30 – 34 ¾	*2
88 – 101	34 ¾ – 39 ¾	*3
101 – 115	39 ¾ – 45 ½	*4
115 – 130	45 ½ – 51 ¼	*5

WHY CHOOSE MANUMED ACTIVE E+MOTION?

Breathable, seamless sports knit and ergonomic stabilising rod (adjustable).

Aid number 05.07.02.0127

WHY CHOOSE LUMBAMED E+MOTION?

Breathable sports knit and adjustable pad for massage and stabilisation.

Aid number 23.14.03.2018





**PUSH
YOUR
PERFOR
MANCE**



PERFORMANCE

TESTED BY

TOP ATHLETES



*From the German Ski Association (DSV) to the Bavarian Olympic Training Centre – top athletes have put the E*motion products through their paces. For more information, simply scan the QR code.*

medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
Phone: +49 921 912-333
Fax: +49 921 912-57
medi@medi.de
www.medi.de

Official supplier of the
German national ski teams



Olympiastützpunkt
Bayern



4 099787 244147

0097EN6 / 03.2023